



promed®

Inspiratory Exerciser

Enabling You To Deliver Quality Care

The promed® 3-ball, flow-oriented incentive spirometer features three colour-coded balls in chambers to measure inspired volume.

The respirometer is used to help improve lung performance in conditions where lung function may be compromised due to various respiratory disorders. This device enables patients to exercise their lungs with the ability to visually examine progress.

promed® Inspiratory Exerciser Features:

- Wide flow rate range from 600 - 1200mL/sec
- Three colour-coded balls/three chambers
- Minimum flow imprinted on each chamber
- Compact design and break-resistant plastic
- Light-weight
- Compact
- Break-resistant plastic
- Removable tubing for cleaning convenience
- Clear tubing enables clinical inspection of any secretions
- Anatomically designed mouthpiece for patient comfort



Code	Product	Packaging
PMD00120	Respirometer	20/pk



ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-292
To fax an order, use 1800-067-639
Visit us online: www.thermofisher.com.au

©2014 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917