




20% off selected Nunc EasYFlask Tissue Culture Flasks

pH testing of the growth media in the culture flasks each day is critical to maintaining the system in a healthy state for cell growth.

The Nunc™ EasYFlask™ with Nunclon™ Delta surface treatment is designed for optimum cell growth. Each batch is tested against 3 cell types and 4 cell lines with Gibco™ media to ensure consistent, reproducible lot to lot quality.

The vent and filter flasks allow researchers to visually check that flasks are set to an open or closed position, and the unique lid requires only a one third turn to allow gas exchange.



Item	Description	Promo Price	Case of #
NUN156340	EasYFlask 25cm ² Nunclon surface with vent/close cap	\$131.65	200
NUN156367	EasYFlask 25cm ² Nunclon surface with filter cap	\$165.00	200
NUN156472	EasYFlask 75cm ² Nunclon surface with vent/close cap	\$125.00	100
NUN156499	EasYFlask 75cm ² Nunclon surface with filter cap	\$138.00	100
NUN159920	EasYFlask 175cm ² Nunclon surface with vent/close cap	\$76.00	30
NUN159910	EasYFlask 175cm ² Nunclon surface with filter cap	\$76.00	30

Find out more at thermofisher.co.nz

*Pricing valid until 18th August 2017. All Prices are in New Zealand dollars and exclude GST. All offers subject to Thermo Fisher standard terms and conditions. Errors and omission exempted. Specifications, terms and pricing are subject to change. Please consult your local sales representative for details. © 2017 Thermo Fisher Scientific Inc. All rights reserved. Trademarks used are owned as indicated on thermofisher.co.nz/trademarks. **1498184162**

Order Placement:

For customer service, call 0800-933-966
 To fax an order, use 0800-329-246
 To email an order, ordersnz@thermofisher.com

Contact Us:

For customer service, call 0800-933-966
 For service and calibration, call 0800-933-966
 To order online: thermofisher.co.nz

ThermoFisher
 SCIENTIFIC