



Azlon laboratory plasticware

Offers on all your plasticware essentials

Baskets

Low form

- Strong lightweight polypropylene rectangular basket with good chemical resistance
- Suitable for use as a draining or autoclave basket and for general storage
- Autoclavable to 121°C
- Able to stack when not in use to save space



Item Code	Description	Unit	Price
AZLTWR200	Basket draining 10 L 400 x 300 x 100 mm (L x W x H)	EA	\$21.90

With hinged lid

- Strong polypropylene basket with hinged lid
- Useful for storage in the laboratory
- Stackable design for space saving when not in use
- Autoclavable at 121°C
- Ideal for sterilising test tubes



Item Code	Description	Unit	Price
AZLBWK006	Basket white PP 167 x 167 x 155 mm (L x W x H)	EA	\$22.90
AZLBWK009	Basket white PP 225 x 225 x 230 mm (L x W x H)	EA	\$24.85

Beakers

Straight sided with printed graduations

- Polypropylene beaker with straight sides and stable shape
- Blue graduations are clearly printed in indelible solvent resistant ink
- Autoclavable at 121°C with good chemical resistance
- Ideal for general laboratory use



Item Code	Description	Unit	Price
AZLBDA218P	Beaker 50 mL PP straight side printed graduations	10 pcs	\$27.30
AZLBDA224P	Beaker 100 mL PP straight side printed graduations	10 pcs	\$31.00
AZLBDA230P	Beaker 250 mL PP straight side printed graduations	10 pcs	\$43.55
AZLBDA236P	Beaker 500 mL PP straight side printed graduations	10 pcs	\$58.65
AZLBDA240P	Beaker 1 L PP straight side printed graduations	5 pcs	\$34.15
AZLBDA242P	Beaker 2 L PP straight side printed graduations	EA	\$12.75

Tapered with printed graduations

- High clarity polypropylene beakers with tapered shape
- Blue, easy to read printed graduations
- Autoclavable at 121°C with good chemical resistance
- Precise non-drip pouring spouts
- Easy to handle and stackable for space saving



Item Code	Description	Unit	Price
AZLBPB0050P	Beaker 50 mL PP tapered blue printed graduations	10 pcs	\$27.10
AZLBPB0100P	Beaker 100 mL PP tapered blue printed graduations	10 pcs	\$32.95
AZLBPB0500P	Beaker 500 mL PP tapered blue printed graduations	10 pcs	\$43.55
AZLBPB1000P	Beaker 1 L PP tapered blue printed graduations	5 pcs	\$25.75
AZLBPB2000P	Beaker 2 L PP tapered blue printed graduations	EA	\$8.55
AZLBPB3000P	Beaker 3 L PP tapered blue printed graduations	EA	\$36.95

Bottles

Narrow neck polypropylene

- Translucent, rigid polypropylene bottles with polypropylene screw caps
- Leakproof, linerless design
- Autoclavable at 121°C with good chemical resistance
- Suitable for food contact applications, liquid storage, media preparation and sampling applications
- Heavy duty and very durable in use



Item Code	Description	Unit	Price
AZLBNP0100PN	Bottle 100 mL PP round narrow neck PP screw cap	10 pcs	\$17.65
AZLBNP0150PN	Bottle 150 mL PP round narrow neck PP screw cap	10 pcs	\$19.45
AZLBNP0250P	Bottle 250 mL PP round narrow neck PP screw cap	10 pcs	\$29.85
AZLBNP0500P	Bottle 500 mL PP round narrow neck PP screw cap	10 pcs	\$36.85
AZLBNP1000P	Bottle 1 L PP round narrow neck PP screw cap	5 pcs	\$28.30
AZLBNP2000P	Bottle 2 L PP round narrow neck PP screw cap	5 pcs	\$61.45

Wide neck polypropylene

- Translucent, rigid polypropylene bottles with polypropylene screw caps
- Smooth inner contours for efficient pouring and ease of cleaning
- Autoclavable at 121°C with good chemical resistance
- Suitable for food contact applications, liquid storage, media preparation and sampling applications
- Leakproof, heavy duty and very durable in use



Item Code	Description	Unit	Price
AZLBLP0100P	Bottle 100 mL PP round wide neck PP cap	10 pcs	\$14.60
AZLBLP0150P	Bottle 150 mL PP round wide neck PP cap	10 pcs	\$16.00
AZLBWP0250P	Bottle 250 mL PP round wide neck PP cap	10 pcs	\$29.85
AZLBWP0500P	Bottle 500 mL PP round wide neck PP cap	10 pcs	\$35.90
AZLBWP1000P	Bottle 1 L PP round wide neck PP cap	5 pcs	\$38.25
AZLBWP2000P	Bottle 2 L PP round wide neck PP cap	5 pcs	\$46.80

Amber wide neck polypropylene

- Opaque, rigid amber polypropylene bottles with black polypropylene screw caps
- Ideal for a wide range of light sensitive liquid or solid storage
- Autoclavable at 121°C
- Leakproof and heavy duty
- Good general chemical resistance



Item Code	Description	Unit	Price
AZLBWP0060AP	Bottle 60 mL PP amber round wide neck PP cap	10 pcs	\$14.80
AZLBWP0150AP	Bottle 150 mL PP amber round wide neck PP cap	10 pcs	\$20.75
AZLBWP0250AP	Bottle 250 mL PP amber round wide neck PP cap	10 pcs	\$32.95
AZLBWP0500AP	Bottle 500 mL PP amber round wide neck PP cap	10 pcs	\$39.55
AZLBWP1000AP	Bottle 1 L PP amber round wide neck PP cap	5 pcs	\$42.30
AZLBWP2000AP	Bottle 2 L PP amber round wide neck PP cap	5 pcs	\$51.75

Square high density polyethylene

- Square high density polyethylene bottles with a plug seal and screw cap
- Heavy duty with an easy to grip cap
- Suitable for food contact and liquid or solid storage
- Moulded graduations for approximate volume measurement



Item Code	Description	Unit	Price
AZLBGE224P	Bottle 100 mL HPDE square wide mouth graduated	10 pcs	\$26.20
AZLBGE230P	Bottle 250 mL HPDE square wide mouth graduated	10 pcs	\$36.25
AZLBGE336P	Bottle 500 mL HPDE square wide mouth graduated	10 pcs	\$45.55
AZLBGE340P	Bottle 1 L HPDE square wide mouth graduated	5 pcs	\$34.80
AZLBGE342P	Bottle 2 L HPDE square wide mouth graduated	5 pcs	\$54.50

Dispensing

- Round low density polyethylene bottles fitted with a spouted dropper closure
- Closure features a captive sealing cap which can be cut to adjust flow rate
- Ideal for use as a dropping bottle and for dispensing liquids
- Good chemical resistance



Item Code	Description	Unit	Price
AZLBGE814PNC	Bottle dispensing 30 mL LDPE spouted cap	10 pcs	\$19.55
AZLBGE820PNC	Bottle dispensing 60 mL LDPE spouted cap	10 pcs	\$22.25
AZLBGE828PNC	Bottle dispensing 150 mL LDPE spouted cap	10 pcs	\$27.95
AZLBGE836PNC	Bottle dispensing 500 mL LDPE spouted cap	10 pcs	\$56.90

Wash bottles

Classic style

- Traditional and square shoulder narrow neck round wash bottles
- Translucent and flexible low density polyethylene with good chemical resistance
- Fitted with one-piece polypropylene closures that maintain their shape in use
- Closures feature precision jet-tip aperture for fine control



Item Code	Description	Unit	Price
AZLWGF230P	Wash bottle 250 mL round LDPE white closure	5 pcs	\$27.90
AZLWGF236P	Wash bottle 500 mL round LDPE white closure	5 pcs	\$33.90

Integral spout

- Slim narrow neck low density polyethylene wash bottles with good chemical resistance
- Ergonomic wash bottles with an integrally moulded dispensing tube jet tip
- Precision tip is ideal for fine liquid jet dispensing



Item Code	Description	Unit	Price
AZLBGF130P	Wash bottle 250 mL integral narrow neck LDPE	EA	\$12.30
AZLBGF136P	Wash bottle 500 mL integral narrow neck LDPE	EA	\$13.00

Multi-lingual GHS labelled

- Wide neck low density polyethylene wash bottles with solvent labels
- Printed in 4 languages (English, French, German and Spanish)
- Solvent labels include GHS symbols, NFPA diamond, risk phrases and safety information
- Fitted with one-piece polypropylene closures that maintain their shape in use
- Closures feature precision jet-tip aperture for fine control



Item Code	Description	Unit	Price
AZLWGW536PML	Wash bottle 500 mL non-venting white water	5 pcs	\$47.95
AZLWGW601PML	Wash bottle 500 mL non-venting white distilled water	5 pcs	\$47.95
AZLWGW603PML	Wash bottle 500 mL non-venting white deionised water	5 pcs	\$47.95
AZLWGW537PML	Wash bottle 500 mL non-venting red acetone	5 pcs	\$47.95
AZLWGW538PML	Wash bottle 500 mL non-venting green methanol	5 pcs	\$47.95
AZLWGW539PML	Wash bottle 500 mL non-venting blue isopropanol	5 pcs	\$47.95
AZLWGW541PML	Wash bottle 500 mL non-venting orange ethanol	5 pcs	\$47.95

Do-it-yourself safety venting

- Do-it-yourself DripLok® wash bottles with blank write-on-panel
- Printed in 3 languages (English, French and German)
- Supplied with adhesive hazard labels



Item Code	Description	Unit	Price
AZLWGW001VTML	Wash bottle 250 mL safety venting white DIY label	5 pcs	\$58.90
AZLWGW002VTML	Wash bottle 500 mL safety venting white DIY label	5 pcs	\$61.00

Connectors

Straight pattern

- White polypropylene straight pattern tubing connectors
- For connecting tubing of similar or different sizes
- Tapered arms are rigid and stepped for a firm grip
- Autoclavable at 121°C with good chemical resistance



Item Code	Description	Unit	Price
AZLCWH600	Connector straight fits 4-7 mm tubing	10 pcs	\$37.25
AZLCWH602	Connector straight fits 4-7 mm and 7-12 mm tubing	10 pcs	\$38.85
AZLCWH604	Connector straight fits 7-12 mm tubing	10 pcs	\$40.90
AZLCWH606	Connector straight fits 4-7 mm and 10-15 mm tubing	10 pcs	\$42.40
AZLCWH608	Connector straight fits 7-12 mm and 10-15 mm tubing	10 pcs	\$44.15
AZLCWH610	Connector straight fits 10-15 mm tubing	10 pcs	\$46.80

“Y” pattern

- Polypropylene “Y” pattern tubing connectors
- Tapered arms are rigid and stepped for a firm grip
- Autoclavable at 121°C with good chemical resistance



Item Code	Description	Unit	Price
AZLCWH304	Connector Y tapered 4-7 mm PP clear	10 pcs	\$46.55
AZLCWH306	Connector Y tapered 6-7 mm PP clear	10 pcs	\$51.20
AZLCWH308	Connector Y tapered 8-9 mm PP clear	10 pcs	\$56.00
AZLCWH310	Connector Y tapered 10-12 mm PP clear	10 pcs	\$59.80

Containers

- Wide neck polypropylene containers with linerless screw caps
- Ideal for storage of powders and samples
- Autoclavable at 121°C with good chemical resistance



Item Code	Description	Unit	Price
AZLJWP0060P	Container 60 mL PP wide neck 53 x 47 mm (OD x height)	48 pcs	\$91.10
AZLJWP0125P	Container 125 mL PP wide neck 70 x 60 mm (OD x height)	36 pcs	\$94.50
AZLJWP0250P	Container 250 mL PP wide neck 89 x 66 mm (OD x height)	36 pcs	\$186.60
AZLJWP0500P	Container 500 mL PP wide neck 89 x 97 mm (OD x height)	24 pcs	\$134.40

Cylinders

Moulded graduations

- Polypropylene measuring cylinders with moulded graduations
- Calibrated to Class B tolerances
- Large pouring spout for easy dispensing
- Stable hexagonal base with anti-suction design that prevents base from sticking to a wet surface
- Good chemical resistance



Item Code	Description	Unit	Price
AZLCP0025P	Cylinder measuring 25 mL PP moulded graduations	10 pcs	\$70.00
AZLCP0050P	Cylinder measuring 50 mL PP moulded graduations	10 pcs	\$80.40
AZLCP0100P	Cylinder measuring 100 mL PP moulded graduations	5 pcs	\$53.25
AZLCP0500P	Cylinder measuring 500 mL PP moulded graduations	EA	\$23.25
AZLCP1000P	Cylinder measuring 1 L PP moulded graduations	EA	\$38.95
AZLCP2000P	Cylinder measuring 2 L PP moulded graduations	EA	\$85.40

Printed graduations

- Polypropylene measuring cylinders with blue printed graduations
- Calibrated to Class B tolerances
- Clear pin sharp blue graduations
- Stable hexagonal base with anti-suction design that prevents base from sticking to a wet surface
- Good chemical resistance



Item Code	Description	Unit	Price
AZLCPB0025P	Cylinder measuring 25 mL PP printed graduations	10 pcs	\$93.90
AZLCPB0050P	Cylinder measuring 50 mL PP printed graduations	10 pcs	\$107.65
AZLCPB0100P	Cylinder measuring 100 mL PP printed graduations	5 pcs	\$64.45
AZLCPB0500P	Cylinder measuring 500 mL PP printed graduations	EA	\$24.40
AZLCPB1000P	Cylinder measuring 1 L PP printed graduations	EA	\$38.25
AZLCPB2000P	Cylinder measuring 2 L PP printed graduations	EA	\$72.15

Funnels

Lightweight

- Lightweight polypropylene funnels with ribs that prevent air locks, allowing for rapid filtration
- Good chemical resistance
- Autoclavable to 121°C



Item Code	Description	Unit	Price
AZLFWC058	Funnel 50 mm diameter PP lightweight 78 mm height	10 pcs	\$28.00
AZLFWC060	Funnel 75 mm diameter PP lightweight 107 mm height	10 pcs	\$25.90
AZLFWC062	Funnel 100 mm diameter PP lightweight 149 mm height	10 pcs	\$53.55
AZLFWC064	Funnel 120 mm diameter PP lightweight 170 mm height	5 pcs	\$35.45
AZLFWC066	Funnel 150 mm diameter PP lightweight 210 mm height	5 pcs	\$45.95

Powder

- Polypropylene powder funnel that are unbreakable in normal use
- Parallel stem minimises 'bridging' of powder and improves flow
- Autoclavable at 121°C



Item Code	Description	Unit	Price
AZLFWC204	Funnel 80 mm diameter PP powder 74 mm height	5 pcs	\$38.10
AZLFWC206	Funnel 100 mm diameter PP powder 90 mm height	5 pcs	\$42.10
AZLFWC208	Funnel 120 mm diameter PP powder 106 mm height	5 pcs	\$49.00
AZLFWC210	Funnel 150 mm diameter PP powder 135 mm height	5 pcs	\$76.80

Ice buckets

- Durable, rigid polyurethane skinned foam ice bucket with lid
- Excellent insulating properties make it suitable for use with water-ice, dry-ice and salt solutions
- Integral handles allow for easy lifting and carrying
- Capacity of 4.5 litres
- Dimensions 185 x 330 mm (H x W) and weight 1200 g



Item Code	Description	Unit	Price
AZLIBB001P	Ice bucket 4.5 L polyurethane 185 x 330 mm black	EA	\$190.00
AZLIBB002P	Ice bucket 4.5 L polyurethane 185 x 330 mm blue	EA	\$190.00
AZLIBB003P	Ice bucket 4.5 L polyurethane 185 x 330 mm green	EA	\$190.00

Jugs

- High clarity polypropylene tapered jugs with moulded graduations
- Suitable for food contact applications
- Autoclavable at 121°C with good chemical resistance
- Permanent moulded graduations that won't rub off
- Ergonomically designed handle



Item Code	Description	Unit	Price
AZLJPM0500P	Jug 500 mL PP with handle moulded graduations	EA	\$17.50
AZLJPM1000P	Jug 1 L PP with handle moulded graduations	EA	\$21.00
AZLJPM2000P	Jug 2 L PP with handle moulded graduations	EA	\$28.05
AZLJPM5000P	Jug 5 L PP with handle moulded graduations	EA	\$37.40

Pipette fillers

Bibbette

- Bibbette® pipette filler makes liquid handling simple and easy and prevents contamination through a replaceable PTFE membrane
- PTFE and silicone rubber valve assembly controls the vacuum created by simply squeezing the unique rubber bulb
- Offers greater control and is able to easily and accurately fill or empty volumes of 0.1 mL to 100 mL
- Blow out pipettes are catered for by the depression of a small thumb operated squeeze ball
- A replaceable PTFE membrane filter ensures liquid is not accidentally drawn into the unit (2 supplied)
- The filter and pipette holding cone can be sterilised by autoclaving



Item Code	Description	Unit	Price
AZLPFL001B	Pipette filler bibbette blue	EA	\$187.40

Natural rubber

- Natural rubber pipette filler for basic liquid handling
- Suitable for pipettes up to 100 mL
- Pipette filler works via a series of glass ball valves for air release, suction and delivery



Item Code	Description	Unit	Price
AZLPWY144	Pipette filler rubber 3 valve standard	EA	\$24.95
AZLPWY145	Pipette filler rubber 3 valve for single handed operation	EA	\$25.30

Draining rack

- Free standing draining rack made from chemically resistant polyamide epoxy coated steel
- Suitable for use with plastic and glass labware
- Coating offers greater heat and chemical resistance than other coatings
- Has eight 110 mm and twenty four 80 mm long pegs. Peg diameter is 4.5 mm
- Overall dimensions of the rack are 410 x 300 mm



Item Code	Description	Unit	Price
AZLRWK070FS	Rack draining with integral free standing support	EA	\$146.20

Stirrer bars with pivot ring

- Strong alloy magnets sheathed in pure white inert polytetrafluoroethylene
- Central pivot ring ensures efficient spinning even on curved or uneven bases
- Excellent chemical resistance and can withstand temperatures up to 260°C



Item Code	Description	Unit	Price
AZLSWN502	Stirrer bar PTFE length 12 mm OD 4.5 mm pivot ring	2 pcs	\$9.00
AZLSWN503	Stirrer bar PTFE length 20 mm OD 6 mm pivot ring	2 pcs	\$11.70
AZLSWN505	Stirrer bar PTFE length 35 mm OD 6 mm pivot ring	2 pcs	\$15.80
AZLSWN507	Stirrer bar PTFE length 40 mm OD 8 mm pivot ring	2 pcs	\$18.45
AZLSWN508	Stirrer bar PTFE length 60 mm OD 9.5 mm pivot ring	2 pcs	\$31.10
AZLSWN510	Stirrer bar PTFE length 70 mm OD 9.5 mm pivot ring	2 pcs	\$29.00

Water jet pump

- Suitable for use with most laboratory water pressures
- Integrated non-return valve increases safety
- Detachable vacuum connection that accepts 6-9 mm bore tubing
- Water inlet hose connection accepts 10-12 mm bore tubing
- Very low water consumption: 190 litre/h at 3.5 bar water supply pressure
- Able to achieve pressure of 16 mbar (Water temperature of 12°C) with a water supply pressure from 3-6 bar
- High suction capacity with a flow rate of approximately 400 litres/h of air (At 1 atm and water supply pressure of 3.5 bar at 12°C)
- Supplied with water inlet hose connection and R 3/4 inch threaded tap connectors, with R 1/2 inch reduction adapter



Item Code	Description	Unit	Price
AZLPWD003	Pump water jet PP with integral non return valve	EA	\$161.95



Order Placement:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

Contact Us:

For customer service, call 0800-933-966
For service and calibration, call 0800-933-966
To order online: thermofisher.co.nz

Find out more at thermofisher.co.nz

ThermoFisher
SCIENTIFIC