



Reduce specimen identification errors Enhance lean workflow

Thermo Scientific™ SlideMate™ AS On-Demand Slide Printer

The Thermo Scientific™ SlideMate™ AS is the next generation Thermo Scientific Slide Printer delivering a smaller footprint than previous models, making it even easier to fit next to the microtome. The SlideMate AS promotes lean workflow through unique features designed specifically for on-demand printing.

- Fully integrated colour touch-screen and simple software are easy to use
- Built-in barcode scanner reads cassette information and imports data to the slide label
- Connects easily to Laboratory Information Systems
- Prints directly onto the slide, eliminating handwriting or labels



Printing using the SlideMate AS is as easy as 1,2,3...

1. Scan the cassette
2. Insert the slide (or use the slide delivery system)
3. Retrieve the printed slide

Talk to your local product specialist about how the SlideMate AS Slide Printer can improve your laboratory workflow today.

Authorised distributor

In Australia:

For customer service, call 1300-735-292
To fax an order, use 1800-067-639
To email an order, ordersau@thermofisher.com

In New Zealand:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

ThermoFisher
SCIENTIFIC