



Science.
Applied to Life.™

3M Food Safety

Bioluminescence. Applied to guarding the line.

3M™ Clean-Trace™ Hygiene Monitoring
and Management System





Ready when you are.

3M™ Clean-Trace™ ATP Tests

You don't have time for inefficient testing — you need to get in, get your sample and get out. 3M Clean-Trace ATP Tests feature an intuitive design that's easy to use and implement, with a long, flexible shape, allowing you to swab even difficult-to-reach areas. And the sleek, portable swabs are easy to activate with one hand.

Flexible design helps you swab difficult areas

Easy to activate for increased efficiency

Pre-moistened and ready to use — just grab and go



Hygiene monitoring has a new hero.



3M™ Clean-Trace™ Luminometer

The newly redesigned 3M Clean-Trace Luminometer features improved ease of use, faster time to result and simple, one-handed operation. Its ergonomic design makes testing simple, minimizing training time and costs. And with a higher degree of repeatability and lower variability, you'll get accurate results you can count on.

Simple, one-handed operation allows multitasking

Wireless capabilities with Bluetooth® Technology

New ergonomic design is comfortable to hold



Powerful. Flexible. Indispensable.

3M™ Clean-Trace™ Hygiene Management Software

Don't waste time searching for reports. Featuring a redesigned dashboard for a quick and easy way to create, store and retrieve reports, 3M Clean-Trace Hygiene Management Software turns data into information, helping you identify root causes and take action quickly. Now even easier to use and more intuitive with convenient wireless and networking capabilities. Put your data to work identifying workplace trends, optimizing productivity and adding value to the process.



Data that works for you.

Data is just data — until it's simplified, streamlined and presented as information that helps you make better, more educated decisions. With 3M Clean-Trace Hygiene Management Software, you can quickly pinpoint exactly where a sanitation problem occurred and address it. Find areas of concern or predict where issues could happen. Confidently reduce unnecessary recleaning and downtime for faster production. Be proactive and prevent small problems before they become big ones.

- ▶ Intuitive, user-friendly software optimizes productivity
- ▶ Redesigned dashboard consolidates information
- ▶ Wireless data transfer streamlines workload
- ▶ Automates data transfer, report generation and delivery

The test that sees beyond the surface.

When people are relying on you to identify and manage risk, you need accurate information. And to confidently make these high-risk decisions, you need more than a visual assessment. You need faster results than microbiological tests can deliver. You need a reliable, complete system to help you manage your hygiene monitoring — so you can make the call and get things moving.

	Visual Assessment	Microbiological Tests	ATP Tests
Rapid	+	-	+
Sensitive	-	+	+
Quantitative	-	+	+
Detect product residues	+	-	+
Simple	+	-	+

Microbes and product residue contain Adenosine Triphosphate (ATP), an indicator of biological residues that can be easily detected to measure cleanliness quickly. And unlike other methods, ATP testing is rapid, sensitive, quantitative, effective and simple.



Support you can count on.

Protecting the world's food supply is a big responsibility — but you don't have to do it alone. We have people to support you locally by providing technical help and exceptional service. It's our job to help you mitigate risk and improve operations — every batch, every shift, every day.

To learn more about how the 3M™ Clean-Trace™ Hygiene Monitoring and Management System can help protect your brand and your business, contact your local 3M representative or visit 3M.com/foodsafety/Clean-Trace.

For more information on the 3M™ Clean-Trace™ Hygiene Monitoring and Management System:

go.3M.com/clean-trace-hygiene-monitoring3



Authorised distributor

Order Placement:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

Contact Us:

For customer service, call 0800-933-966
For service and calibration, call 0800-933-966
To order online: thermofisher.co.nz

Ordering Information

▶ 3M Clean-Trace Luminometer and Data Trending Software*

Product #	Description
LM1	1/each

*Computer not included.

▶ 3M Clean-Trace Surface ATP Test


Product #	Description
UXL100	100/case

▶ 3M Clean-Trace Water Plus — Total ATP Test

Product #	Description
AQT200	100/case

▶ 3M Clean-Trace Water-Free ATP Test

Product #	Description
AQF100	100/case



ThermoFisher
SCIENTIFIC

© 3M 2016. All rights reserved. 3M and Clean-Trace are trademarks of 3M. All other trademarks are the property of their respective owners. Used under license in Canada. Please recycle. Printed in U.S.A.

70-2011-5084-7